


H/I Cryo-Histomat® MK-2



H/I Cryo-Histomat MK2 shown with H/I Sliding Microtome

THERMO-ELECTRONIC MICROTOME FREEZING ATTACHMENT

Features and Benefits:

- Converts any standard microtome for frozen sectioning
- "Solid State" - no moving parts, no wear, high reliability
- Fast freezing preserves cell detail
- No circulating refrigerants or cumbersome CO₂ cylinders
- Positive thaw control for quick specimen change
- Section temperatures to -40° C
- Low running cost - about 2 cents per hour
- Noiseless - vibration free
- Compact - conveniently portable

***CRYO-HISTOMAT MK-2** is a "new technology" freezing attachment for any microtome. Fundamentally different from other freezers, **Cryo-Histomat** dispenses with bulky, inconvenient refrigerants. Freezing is accomplished and maintained by the "Peltier Effect" in a compact, easy-to-operate thermo-electric device.*

Installation

Quick and easy ... just place the cryo-stage in the clamp of your microtome, plug into mains outlet and connect the plastic tube to a lab water faucet. Alternately, in areas remote from water supply, use the H/I closed circuit pump and reservoir. Only a moderate flow of water is needed to carry off heat removed from the specimen by the Peltier module.

Operation

Simple and convenient ... Place your specimen on the Peltier stage suitably orientated for the microtome, switch on and adjust the current control until meter reads 10 amps. Freezing at the specimen/stage interface is complete within seconds. If reorienting is necessary, or to remove tissue block after sectioning is completed, depress the spring toggle switch to "thaw" and hold. The block will free and refreeze to the stage when the switch is released. Freezing can be maintained indefinitely.

Safety is Built In

Only low voltage is transmitted to the cryo-stage; it never exceeds 25 v. The stage is protected from overheating by thermal switch within the module which cuts off power to the solid-state devices well before a destructive temperature is reached.



SPECIFICATIONS

CONTROL UNIT

Compact, contemporary design. Metal Chassis and case with two carrying handles and rubber feet. Finished in charcoal gray textured vinyl with light blue facia. On-off and freeze-thaw switches, red indicator light, current-temperature control and ammeter 0 - 15. Six (6) foot mains supply cord with molded 3-pin plug.

Input: 110-125 or 220v., A.C. Output to stage: Maximum 25v., D.C

OPERATING INSTRUCTIONS

Supplied for routine histology - plus suggestions for best temperatures and advanced applications.

CRYO-STAGE

"Solid State" exchanger employing "Peltier Effect"; special design maintains integrity under conditions of extreme thermal expansion/contraction; DC low voltage connector; soft vinyl tubing for water supply.

Stage Area Model A: 32 x 32mm - larger stage area available for sliding microtomes (see ordering information below)

FREEZING RANGE

To -40° C. under normal ambience of water and laboratory environment

MOUNTING HARDWARE

Three adaptors; block 19mm square; two short pillars 9.5 and 120mm diameter; for all standard microtomes.

PHYSICAL DATA

Width: 33 cm (13") plus 2 handles 6 cm; Depth: 25.5 cm (10"); Height: 16.5 cm (6.5"); Weight: 9 kilos net, 12 kilos gross.

ACKNOWLEDGEMENT: to Jean Peltier, who in 1834 discovered the principle of thermoelectric cooling. Utilization became possible only with spin-off from transistor technology.

Ordering Information:

MANUAL TEMPERATURE CONTROL

Models A and B control units are identical in appearance and operating Model B has increased output to operate larger stage.

- | | Cat.No. |
|---|----------------|
| • Model A, 30 x 40 mm cryo-stage - recommended for clinical use on standard rotary, sliding and freezing microtomes. 110-125v 50/60Hz | 11 00 01 |
| • Model B, 75 x 80 mm cryo-stage - for larger sliding microtomes | 11 00 02 |

AUTOMATIC TEMPERATURE CONTROL

Generally as manual control Cryo-Histomat but with specimen temperature feed-back via microprobe needle for automatic control. Desired temperature may be set in the range -40° C and read to 0.1° C. Probe and readout unit may be used independently as laboratory thermometer. Probe supplied 24 or 30 gauge. 110-125v; 50/60Hz

- | | |
|---|----------|
| • Model A auto, 30 x 40 mm cryo-stage | 11 00 03 |
| • Model B auto, 75 x 80 mm cryo-stage | 11 00 04 |

MICROPROBE THERMOMETER

Use with manual control Cryo-Histomat to show temperature of tissue specimen - or as laboratory thermometer Microprobe thermometer; temperature range - 100° to + 170° in 5 steps; reads to 0.1° C.

Probe supplied 23 gauge, 3cm long

11 00 12

CLOSED CIRCUIT COOLANT SYSTEM

Allows Cryo-Histomat operation remote from water supply or where water must be conserved.

Water pump and 2 liter reservoir 110-125v, 60 cycle operation

11 00 10

HACKER Instruments & Industries Inc

1132 Kincaid Bridge Road • P.O. Box 1176

Winnsboro, South Carolina 29180 USA

Tel: (803) 712-6100 • Fax: (803) 712-6116

www.HackerInstruments.com • HACKERLAB@aol.com