

Bright

2212

from

Hacker

The Only Truly Portable Cryostat

*Ideal for labs with limited space or
as a field cryostat*

Complete with: -30° operation, digital temperature displays for chamber and freezer; auto-defrost with buffered quartz clock; special microtome; 2/12 µm advance; anti-roll device; magnetic blade holder; 2 blades; 2 object holders; spare fuses



*Dimensions: 28.5" (L) x 17.2" (W) x 11.8" (H)
(740 mm x 545 x 315)*

*Shipping Weight: 90 lbs (37 kg)
110-125V, 60 Hz, hospital grade plug
(220V also available)*

Hacker-Bright MicroCryostat 2212

Cat. No. 11 50 04

	Cat. No.	Ancillaries
	11 55 06	Quick-freeze attachment (CO ₂ required)
	11 55 12	Cryo-mold, 22 mm diameter
	11 55 14	Cryo-mold, 16 mm diameter
	11 55 16	Cryo-mold, 10 mm diameter
<i>Options</i>	11 55 22	U.V. Sterilization unit
	11 55 24	Socket for orientable object holders
	11 55 26	Object holder, 22 mm diameter, orientable
	11 55 23	Storage tray for 10 object holders
<i>Extras & Replacements</i>	11 55 04	Object holder, 22 mm diameter
	11 55 18	Anti-roll plate, replacement for standard unit
	11 55 10	Blades, packet of 10**
	10 70 62	"Magnaplate" anti-roll plate
	11 70 76	Cryo-M-Bed, 113 ml (4oz)

Please call 1-800-4-HACKER for current pricing

Hacker
Instruments & Industries Inc.
Since 1942

www.hackerinstruments.com
www.tidylabstore.com

1132 Kincaid Bridge Road • PO BOX 1176
Winnsboro, SC 29180 • (803) 712-6100