


AutoSlipper™

Superior Solutions for CoverSlipping



AUTO SLIPPER™ MODEL# 10 75 70 SHOWN WITH MICROSLIP™ DISPOSABLE COVER GLASS CARTRIDGES & BATTERIES

For Spotless Coverslipping The **HACKER-AutoSlipper™** Automates Hand Coverslipping, Anytime - Anyplace

FAST ★ EASY ★ PORTABLE ★ ECONOMICAL

The **AutoSlipper™** is a revolutionary, *patented* device that automatically dispenses glass cover slips, one at a time, greatly facilitating and expediting the hand coverslipping process.

The **AutoSlipper™** eliminates finger prints on the cover slip and the sticking together of cover glasses. The cover slip "*waste factor*" is reduced to zero, as it eliminates the need of hand coverslipping technicians to retrieve cover slips from conventional packaging, causing breakage and exposure to xylene/solvents.

The **AutoSlipper™** uses disposable cartridges of **MicroSlip™** cover glass, available in all standard sizes. **MicroSlip™** is a pre-cleaned, StickFree cover glass, packaged and sealed in a moisture and debris resistant, disposable cartridge.

The **AutoSlipper™** weighs under 14 ounces fully loaded with cover glass and four AA batteries. The **AutoSlipper™** is truly portable to accommodate multiple coverslipping stations as needed. Ideal for routine coverslipping, frozen sections and special stains.

When ordering, please request: Cat.#

Cat.#

HACKER-AutoSlipper

10 75 70

MicroSlip Disposable Cover Glass Cartridge

10 75 72

Available exclusively through your source for innovation in histopathology:

HI HACKER Instruments & Industries Inc
1132 Kincaid Bridge Road • P.O. Box 1176
Winnsboro, South Carolina 29180, USA
Tel: (803) 712-6100 • Fax: (803) 712-6116
e-mail: HACKERLAB@AOL.COM