

Bright

CliniRF


Ultra Low Rapid Freezer Unit from

Hacker



The CliniRF is a small, free standing, top opening, ultra low temperature freezer unit, designed to rapidly freeze tissue specimens for sectioning or processing. The CliniRF eliminates the need to use CO₂ and LN₂, thus also eliminating the need to purchase and store these substances

- Rapid Freezer Tank to -80° C
- -30° C Freezer area for short term storage of specimens
- Removable block & Object holder adaptor
- Included: 4 oz. Cryo-M-Bed, Cooling block removal tool, Allen Key, spare set of fuses, Operators manual



Freezing Bath
Diameter: 275mm (10.83")
Depth: 136mm (5.35")

Rapid Freezer
Diameter: 95mm (3.75")
Depth: (3.75")

Specifications

Weight: 90 lbs (41 kg); Dimensions: 40" (H) x 18" (W) x 20" (D) (1000 x 450 x 500 mm)

Control Panel



CliniRF Ultra Low Rapid Freezer Unit	Cat. No. 18 00 01
CliniRF Optional Rapid Freeze Timer (7-day Programmable)	Cat. No. 18 00 05
CliniRF Optional Variable Temperature (above -80° C)	Cat. No. 18 00 07
CliniRF Extra-Low Temperature (-100° C)	Cat. No. 18 00 03

Hacker
Instruments & Industries Inc.
Since 1942

www.hackerinstruments.com
www.tidylabstore.com

1132 Kincaid Bridge Road • PO BOX 1176
Winnsboro, SC 29180 • (803) 712-6100