

H/I HACKER PD-2.0 PARAFFIN DISPENSER



HIGH PARAFFIN FLOW RATE

**FOR EMBEDDING LARGER
PARAFFIN BLOCKS...**

*THE
HACKER PD-2.0
Paraffin Dispenser*

**YOUR DEPENDABLE,
USER FRIENDLY,
SECONDARY PARAFFIN SOURCE**

- ◆ 1.75 GALLON CAPACITY
- ◆ STAINLESS STEEL INNER RESERVOIR. RIGID METAL CASTING, HARD COAT FINISH AND DROP ON STAINLESS STEEL LID FOR EASY CLEANING
- ◆ FILTER SCREEN PREVENTS DELIVERY TAP BLOCKAGE
- ◆ DRIP FREE, HIGH FLOW PARAFFIN DELIVERY
- ◆ IMPROVED ALUMINUM IMMERSION HEATER FOR MORE EFFICIENT THERMAL CONDUCTIVITY
- ◆ TWO POSITION LEVER TAP LOCKS IN OPEN FOR TOTAL FLOW CONTROL

Specifications

DIMENSIONS	EXTERNAL DIAMETER : 10" X 10" X 18.5" (25CM X 25CM X 47CM) TAP FROM BENCH: HEIGHT 6" (15CM)	
TEMPERATURE RANGE	AMBIENT TO 80°C (± 2°C)	
POWER SUPPLY/RATING	110 VOLTS 60 HZ 420 WATTS	• 220 VOLTS 50/60 HZ 420 WATTS
CAPACITY	1.75 GALLONS / 6.5 LITERS	
WEIGHT	18 LBS / 8.3 KG (110 VERSION)	• 26 LBS / 11.8 KG (220 VERSION)

Ordering Information HACKER PD-2.0 PARAFFIN DISPENSER 220 VAC, 50/60 Hz 54 00 06
 HACKER PD-2.0 PARAFFIN DISPENSER 110 VAC, 60 Hz 54 00 05

HACKER Instruments & Industries Inc

1132 Kincaid Bridge Road • P.O. Box 1176
 Winnsboro, South Carolina 29180, USA
 Tel: (803) 712-6100 • Fax: (803) 712-6116
 e-mail: HACKERLAB@AOL.COM
<http://www.hackerinstruments.com>